



Diathermy Forceps, Cables & Adaptors



www.elitemedical.com.au

REUSABLE DIATHERMY FORCEPS

This range of reusable Bipolar and Monopolar forceps consists of high quality instruments suitable for a wide variety of surgical procedures. Based on the requirements of leading surgeons and specialists, these forceps are designed and manufactured for optimal clinical performance.

Design & Quality

- Light closing force to reduce hand fatigue.
- Precision tip alignment for improved control.
- High-quality nylon insulation resists abrasion.
- Designed to withstand up to 50 cycles.



Compatibility

- Available to suit all common cables and connections.
- Compatible with all generators manufactured to International Standards (IEC 60601).

Diathermy Cables & Adaptors

The diathermy cables and adaptors range from Elite Medical consists of high quality reusable and single use devices for a variety of surgical procedures. They are designed for use with all common bipolar and monopolar instruments and electro-surgical generators.

Tip Styles



Also available

Some styles and patterns of diathermy forceps are available with Non-Stick tips and have black insulation for easy identification. Non-Stick tips are coated with a unique metal alloy which assists with heat dispersion and reduces adhesion of coagulated tissue.

Custom Made Forceps also available

Forceps shown in this brochure are common sizes and styles used in Australia, and are generally available from stock. However, we are able to custom make diathermy

forceps to suit any particular requirements.

Modifications include:

- Shorter or longer lengths
- Larger tip widths
- Toothed and serrated tips
- Alternative cable connections (i.e. Screw, 4mm Socket, etc)

For further information please talk to one of our customer service team on **1800 354 836**

NOW YOU CAN ORDER ONLINE
www.elitemedical.com.au

Exactly the specialised medical supplies you want, when you need them.

BIPOLAR FORCEPS

Reusable (sterilisable) bipolar forceps. Supplied without cable. Cables available to suit.



**Jewellers Straight Bipolar Forceps
with Stainless Steel Micro Tip**

Code	Length	Tip	Style
ALB-341	115mm	0.5mm	Smooth



**Jewellers Bipolar Forceps
with Stainless Steel Angled Micro Tip**

Code	Length	Tip	Style
ALB-342	115mm	0.5mm	Smooth



**Adson Bipolar Forceps
with Smooth or Serrated Tip**

Code	Length	Tip	Style
ALB-251	120mm	0.5mm	Smooth
ALB-252	120mm	1.0mm	Smooth
ALB-253	120mm	1.5mm	Smooth
ALB-261	120mm	0.5mm	Serrated
ALB-262	120mm	1.0mm	Serrated
ALB-263	120mm	1.5mm	Serrated

BIPOLAR FORCEPS

Reusable (sterilisable) bipolar forceps. Supplied without cable. Cables available to suit.



Straight Bipolar Forceps with Smooth Tip

Code	Length	Tip	Style
ALB-161	150mm	0.5mm	Smooth
ALB-162	150mm	1.0mm	Smooth
ALB-163	150mm	1.5mm	Smooth
ALB-171	175mm	0.5mm	Smooth
ALB-172	175mm	1.0mm	Smooth
ALB-173	175mm	1.5mm	Smooth
ALB-181	200mm	0.5mm	Smooth
ALB-182	200mm	1.0mm	Smooth
ALB-183	200mm	1.5mm	Smooth



Bayonet Bipolar Forceps with Smooth Tip

Code	Length	Tip	Style
ALB-471	175mm	0.5mm	Smooth
ALB-472	175mm	1.0mm	Smooth
ALB-473	175mm	1.5mm	Smooth
ALB-481	200mm	0.5mm	Smooth
ALB-482	200mm	1.0mm	Smooth
ALB-483	200mm	1.5mm	Smooth
ALB-491	225mm	0.5mm	Smooth
ALB-492	225mm	1.0mm	Smooth
ALB-493	225mm	1.5mm	Smooth



Bayonet Bipolar Forceps with Angled Smooth Tip



Code	Length	Tip	Style
ALB-476	175mm	0.5mm	Smooth
ALB-477	175mm	1.0mm	Smooth
ALB-478	175mm	1.5mm	Smooth
ALB-486	200mm	0.5mm	Smooth
ALB-487	200mm	1.0mm	Smooth
ALB-488	200mm	1.5mm	Smooth
ALB-496	225mm	0.5mm	Smooth
ALB-497	225mm	1.0mm	Smooth
ALB-498	225mm	1.5mm	Smooth

MONOPOLAR FORCEPS

Reusable (sterilisable) monopolar forceps. Supplied with a 2.4mm female pin or 4mm banana plug fitting. Supplied without cable. Cables available to suit.





Adson Monopolar Forceps with Serrated Tip

Code 2.4mm Pin	Code 4mm Banana	Length	Tip	Style
				
ALM-261P	ALM-261B	120mm	0.5mm	Serrated
ALM-262P	ALM-262B	120mm	1.0mm	Serrated



Straight Monopolar Forceps with Serrated Tip



Code 2.4mm Pin	Code 4mm Banana	Length	Tip	Style
				
ALM-561P	ALM-561B	150mm	0.5mm	Serrated
ALM-562P	ALM-562B	150mm	1.0mm	Serrated
ALM-563P	ALM-563B	150mm	1.5mm	Serrated
ALM-571P	ALM-571B	175mm	0.5mm	Serrated
ALM-572P	ALM-572B	175mm	1.0mm	Serrated
ALM-573P	ALM-573B	175mm	1.5mm	Serrated
ALM-581P	ALM-581B	200mm	0.5mm	Serrated
ALM-582P	ALM-582B	200mm	1.0mm	Serrated
ALM-583P	ALM-583B	200mm	1.5mm	Serrated

MONOPOLAR FORCEPS

Reusable (sterilisable) monopolar forceps. Supplied with a 2.4mm female pin or 4mm banana plug fitting. Supplied without cable. Cables available to suit.





Bayonet Monopolar Forceps with Serrated Tip

Code 2.4mm Pin	Code 4mm Banana	Length	Tip	Style
				
ALM-481P	ALM-481B	200mm	0.5mm	Serrated
ALM-482P	ALM-482B	200mm	1.0mm	Serrated
ALM-483P	ALM-483B	200mm	1.5mm	Serrated
ALM-491P	ALM-491B	225mm	0.5mm	Serrated
ALM-492P	ALM-492B	225mm	1.0mm	Serrated
ALM-493P	ALM-493B	225mm	1.5mm	Serrated



Bayonet Monopolar Forceps with Angled Serrated Tip

Code 2.4mm Pin	Code 4mm Banana	Length	Tip	Style
				
ALM-486P	ALM-486B	200mm	0.5mm	Serrated
ALM-487P	ALM-487B	200mm	1.0mm	Serrated
ALM-488P	ALM-488B	200mm	1.5mm	Serrated
ALM-496P	ALM-496B	225mm	0.5mm	Serrated
ALM-497P	ALM-497B	225mm	1.0mm	Serrated
ALM-498P	ALM-498B	225mm	1.5mm	Serrated

REUSABLE BIPOLAR CABLES

Universal Bipolar Cables

Eliminate the need to have multiple cables to suit different types of Bipolar Forceps.

- One cable to fit both US (2-pin) Bipolar Forceps and European (Block End) Bipolar Forceps
- Moulded Fixed Pin Generator Plug to meet regulatory requirements
- Available in 4 bright colours for easy identification
- Highly flexible silicone cable
- Reusable and Autoclavable
- Supplied with autoclavable tag and 'pop' punch to track usage.

Code	Colour	Length
V-30/O	Orange	4m
V-30/F	Fuchsia	4m
V-30	Grey	4m
V-30/B	Blue	4m



'Pop' punch system



Available in 4 colours

Bipolar Cables Fly Leads

Code	Colour	Length
S-30	Grey	4m



V-30/B



S-30

REUSABLE BIPOLAR CABLES

Bipolar Cables

Code	Description	Size
190-001-3	Bipolar Cable, Round Angled USA Conn. Shrouded. Fixed Pin	3m
V-16A	Bipolar Cable, Square Angled USA Conn. Not Shrouded. Fixed Pin	3m
V-16/G	Bipolar Cable, Round Straight USA Conn. Not Shrouded. Fixed Pin	3m
S-16A	Bipolar Cable, Square Angled USA Conn. Not Shrouded. Fly Leads	3m
S-16/G	Bipolar Cable, Round Straight USA Conn. Not Shrouded. Fly Leads	3m
V-12	Bipolar Cable, Square Straight EURO Conn. Not Shrouded. Fixed Pin	3m
S-12	Bipolar Cable, Square Straight EURO Conn. Not Shrouded. Fly Leads	3m



V-12



S-12



190-001-3



V-16/G



S-16/G

REUSABLE MONOPOLAR CABLES

Universal Monopolar Cables

Available in 4 colours for easy identification. Fits all common monopolar forceps and comes complete with autoclaveable tag and 'pop' punch to track usage. Double insulated with contrasting inner insulation colour for easy wear/damage identification.

Code	Colour	Size
C-50-O	Orange	8mm plug 4m long
C-50-F	Fuchsia	8mm plug 4m long
C-50-P	Purple	8mm plug 4m long
C-50-B	Blue	8mm plug 4m long



Endoscopic 4.0mm



Foot Controls 4.8mm



Forceps 4.5-5.0mm



UK 4.0mm




Available in 4 colours



'Pop' punch system



C-50-P



REUSABLE MONOPOLAR CABLES

Monopolar Cables

Code	Connection	Length
ALC-623	2.4mm Pin	3m
ALC-624	4mm Banana Plug	3m
ALC-625	4mm Socket	3m
ALC-626	5/32in Screw Fitting	3m
ALC-627	2.5mm Pin. Push/Pull fitting. For snares	3m
C-7/5/R	4mm Shrouded Socket - Red	5m
C-7/5/G	4mm Shrouded Socket - Grey	5m
C-1/R	4.8mm Shrouded Socket - Red	5m



C-1/R



C-7/5/R



ALC-623



ALC-624



ALC-625



ALC-626

SINGLE USE CABLES

Single Use Bipolar Cables

Code	Description	Length
190-014-3	US, 2 Pin, Round, Shrouded	3m
190-014-3ST	US, 2 Pin, Round, Angled	3m
SW46211	US, 2 Pin, Square, Angled	3m

Single Use Monopolar Cables

To suit instruments with 4mm banana fitting

Code	Description	Length
220-024-3	8mm Bovie Plug, Round, Shrouded	3m
220-024-3ES	8mm Bovie Plug, Round	3m
SW47311	8mm Bovie Plug, Round	3m

ADAPTORS

Bipolar Adaptors

Code	Connection
ALC-695	Suits Con-Med style generators
ALC-696	Suits most Erbe style generators

Monopolar Adaptors

Code	Connection
ALC-634	Suits 4mm Banana Pug to 8mm Bovie Plug



**ASK ABOUT OUR
DISPOSABLE
RANGE**



190-014-3



SW46211



220-024-3



SW47311

 *Single use only*



The perfect solution at your fingertips 24/7

Our extensive and trusted range of medical products will fulfil any requirement you may ever have while delivering improved patient outcomes, ease of use and superior value. We carefully source premium instruments from only the most respected manufacturers worldwide. Call our Customer Hotline on 1800 354 836 to see how we can help you.

Electrosurgery

- Electrodes
- Pencils
- Bipolar Forceps, Cables & Adaptors
- Monopolar Forceps, Cables & Adaptors
- Insulated Speculums
- Generators & Smoke Evacuators
- Quivers & Holsters

Surgical Instruments

- Scalpels and Handles
- Scissors
- Forceps
- Suction Tubes
- Retractors & Skin Hooks
- Needle Holders
- Chisels and Gouges
- Fibre Optic Cables
- Laparoscopic Consumables

Instrument Care & Sterilisation

- Instrument Protection Caps
- Identification Tags
- Identification Tape
- Cleaning Brushes
- Cleaning Chemicals & Solutions

Patient Care & Safety

- Pressure Protection
- Operating Table Accessories
- Anaesthetics
- Workplace Health & Safety

About Us

Elite Medical is a family owned and operated company who, since October 1986, have been dedicated to serving the healthcare sector by providing high quality, cost effective and environmentally responsible surgical instruments, equipment and consumables.

Our Values

QUALITY - We proactively seek better results through the continuous improvement of products, processes and people.

WELL-BEING - We promote a safe and supportive atmosphere that provides quality of life for both ourselves and our families.

INTEGRITY - We conduct every transaction with honesty and transparency, standing accountable to our stakeholders and ourselves.

Our People

We believe our people are the most important part of our business. Therefore our team is made up of the most passionate specialists in their field. Our business development managers are committed to searching out innovation and our customer service team go the extra mile to not only meet but exceed your expectations.

Our Community

Elite Medical is 100% Australian owned and operated so we have a practical understanding of your needs. In addition, we are dedicated to being good corporate citizens. In order to achieve this we regularly support charities such as Careflight and the Royal Flying Doctor Service. And we are always improving our practices to ensure they are environmentally responsible.

Exactly the specialised medical supplies
you want, when you need them.



NOW YOU CAN
ORDER ONLINE
www.elitemedical.com.au

Customer Hotline 1800 354 836

P: (07) 3899 1300 **F:** (07) 3399 5858 **E:** sales@elitemedical.com.au

www.elitemedical.com.au