




Endoscope Cleaning Brushes



www.elitemedical.com.au



EBF-1003B Single-Ended Channel Cleaning Brush
Diam: 6mm - Length: 230cm - Box Qty: 50




EBF-1004 Dual-Ended Channel Cleaning Brush
Diam: 6mm - Length: 230cm - Box Qty: 50



EBF-1008 Multipoint Channel Cleaning Brush
Diam: 6mm - Length: 230cm - Box Qty: 50



EBF-1011 Double-ended Combination Valve/Channel Brush
Diam: 11mm/6mm - Length: 230cm - Box Qty: 50



EBF-1015 Double-ended Sweeper Channel Brush
Length: 150mm - Diam: 5mm - Overall Length: 230cm - Box Qty: 50

Endoscopy equipment is sensitive and expensive. Selecting the right endoscope cleaning brush is critical to ensure safe and effective cleaning.

For example, using slightly oversized brushes bend the bristles over, essentially eliminating mechanical abrasion. Brushes with too small a diameter fail to come into firm contact with the cleaning surface, leaving contaminants behind.

Only perfectly sized brushes adequately remove debris and guarantee optimum endoscopy function. Wire ends on channel brushes have plastic tips to prevent endoscope damage.

All brushes are supplied individually packed, non-sterile.


**NOW YOU CAN
ORDER ONLINE**
www.elitemedical.com.au



EBF-1015B Combo Kit, Sweeper Channel Brush. Valve Brush 11/5mm. Length: 150mm - Diam: 5mm - Overall Length: 230cm - Box Qty: 50



EBF-1024 Combo Kit, DE Channel Brush. Diam: 6mm - Length: 230cm. Valve Brush 11/5mm - Box Qty: 50



EBF-1026 Combo Kit, Channel Brush. 1031 Valve Brush 11/5mm. Diam: 3mm - Length: 120cm - Box Qty: 50



EBF-1030 Valve Brush, Double-Ended, Yellow Plastic Handle. Box Qty: 100



EBF-1031 Valve Brush, Double-Ended, Light Blue Plastic Handle. Box Qty: 100



EBF-1032 Valve Brush, Double-Ended, Dark Blue Plastic Handle. Box Qty: 100



EBF-1033 Valve Brush, Single-Ended, Red Plastic Handle. Box Qty: 100



EBF-1034 Valve Brush, Single-Ended, Green Plastic Handle. Box Qty: 100



The perfect solution at your fingertips 24/7

Our extensive and trusted range of medical products will fulfil any requirement you may ever have while delivering improved patient outcomes, ease of use and superior value. We carefully source premium instruments from only the most respected manufacturers worldwide. Call our Customer Hotline on 1800 354 836 to see how we can help you.

Electrosurgery

- Electrodes
- Pencils
- Bipolar Forceps
- Bipolar Cables & Adaptors
- Monopolar Forceps
- Monopolar Cables & Adaptors
- Insulated Speculums
- Generators & Smoke Evacuators
- Quivers & Holsters

Surgical Instruments

- Scalpels and Handles
- Scissors
- Forceps
- Suction Tubes
- Retractors & Skin Hooks
- Needle Holders
- Chisels and Gouges
- Fibre Optic Cables
- Laparoscopic Consumables

Instrument Care & Sterilisation

- Instrument Protection Caps
- Identification Tags
- Identification Tape
- Cleaning Brushes
- Insulation Testers
- Cleaning Chemicals & Solutions

Patient Care & Safety

- Pressure Protection
- Operating Table Accessories
- Anaesthetics
- Workplace Health & Safety

About Us

Elite Medical is a family owned and operated company who, since October 1986, have been dedicated to serving the healthcare sector by providing high quality, cost effective and environmentally responsible surgical instruments, equipment and consumables.

Our Values

INTEGRITY - We act with honesty and transparency, always doing what is right by our stakeholders.

SERVICE - Our customers are our number one priority, therefore we respond to your enquiry with enthusiasm, professionalism and efficiency.

PROGRESS - We are continuously improving our processes in order to remain at the forefront of our industry.

Our People

We believe our people are the most important part of our business. Therefore our team is made up of the most passionate specialists in their field. Our business development managers are committed to searching out innovation and our customer service team go the extra mile to not only meet but exceed your expectations.

Our Community

Elite Medical is 100% Australian owned and operated so we have a practical understanding of your needs. In addition, we are dedicated to being good corporate citizens. In order to achieve this we regularly support charities such as Careflight and the Royal Flying Doctor Service. And we are always improving our practices to ensure they are environmentally responsible.

Exactly the specialised medical supplies
you want, when you need them.



NOW YOU CAN
ORDER ONLINE
www.elitemedical.com.au

Customer Hotline 1800 354 836

P: (07) 3899 1300 **F:** (07) 3399 5858 **E:** sales@elitemedical.com.au

www.elitemedical.com.au

Price List

Endoscope Cleaning Brushes

Version 11014 | Page 1

Order Code	Item Description	Price	Price Each
EBF-1003B	Single-Ended Channel Cleaning Brush, 6mm, 230cm, Box/50	\$95.00	\$1.90
EBF-1004	Dual-Ended Channel Cleaning Brush, 6mm, 230cm, Box/50	\$105.00	\$2.10
EBF-1008	Multipoint Channel Cleaning Brush 6mm, 230cm, Box/50	\$120.00	\$2.50
EBF-1011	Double-ended Combination Valve/Channel Brush, 11mm/6mm, 230cm, Box/50	\$165.00	\$3.30
EBF-1015	Double-ended Sweeper Channel Brush, 150mm/5mm, 230cm, Box/50	\$165.00	\$3.30
EBF-1015B	Combo Kit, Sweeper Channel Brush 150mm/5mm, 230cm, Valve Brush, 11/5mm, Box/50	\$205.00	\$4.10
EBF-1024	Combo Kit, DE Channel Brush 6mm/230cm, Valve Brush 11/5mm, Box/50	\$175.00	\$3.50
EBF-1026	Combo Kit, Channel Brush 3mm/120cm, 1031 Valve Brush 11/5mm, Box/50	\$175.00	\$3.50
EBF-1030	Valve Brush, Double-Ended, Yellow Plastic Handle, Box/100	\$110.00	\$1.10
EBF-1031	Valve Brush, Double-Ended, Light Blue Plastic Handle, Box/100	\$170.00	\$1.70
EBF-1032	Valve Brush, Double-Ended, Dark Blue Plastic Handle, Box/100	\$170.00	\$1.70
EBF-1033	Valve Brush, Single-Ended, Red Plastic Handle, Box/100	\$140.00	\$1.40
EBF-1034	Valve Brush, Single-Ended, Green Plastic Handle, Box/100	\$140.00	\$1.40

Exactly the specialised medical supplies
you want, when you need them.

August 2014. All prices are GST exclusive.



NOW YOU CAN
ORDER ONLINE
www.elitemedical.com.au

Customer Hotline 1800 354 836

P: (07) 3899 1300 F: (07) 3399 5858 E: sales@elitemedical.com.au

www.elitemedical.com.au